



Key Features

- Frequency Conversion Technology
- Maximum Energy Efficiency
- Precise Temperature Control

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l959bp.html

Specifications

Technical Data

Type	-86°C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	HC
Noise Level (dB)	47

Performance

Cooling performance(°C)	-86
Temperature Range(°C)	-40~-86

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
--------------------	-----------------------

Power(W)	1300
Electrical Current(A)	6
Dimensions	
Capacity(L/Cu.Ft)	959/33.9
Net/Gross Weight(approx)	450/485(kg) 992.9/1070.2(lbs)
Interior Dimensions(W*D*H)	1016*716*1310(mm) 40.0*28.2*51.6(in)
Exterior Dimensions(W*D*H)	1296*998*1980(mm) 51.0*39.3*78.0(in)
Packing Dimensions(W*D*H)	1365*1105*2150(mm) 53.7*43.5*84.6(in)
Container load (20'/40'/40'H)	8/16/16
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	Y
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
CO2 Backup System	
LN2 Backup System	
Chart Recorder	

Others

Certificate

CE
