

Reaction Vessels, Cylindrical, Flat Bottom, Graduated

Catalogue Number	Flange Size (DN)	Capacity (ml)	Outer Diameter (mm)	Height (mm)
454 01 01	60	500	65	255
454 01 02	60	1000	105	155
454 02 01	100	500	105	155
454 02 02	100	1000	105	200
454 02 03	100	2000	105	350
454 02 04	100	3000	160	285
454 02 05	100	4000	160	340
454 03 01	150	2000	160	225
454 03 02	150	3000	160	280
454 03 03	150	4000	160	335
454 03 04	150	5000	160	390
454 03 05	150	6000	180	445
454 03 06	150	10000	250	410
454 03 07	150	20000	250	630



Reaction Vessels, Cylindrical, Round Bottom, Graduated

Catalogue Number	Flange Size (DN)	Capacity (ml)	Outer Diameter (mm)	Height (mm)
455 01 01	60	500	65	275
455 01 02	60	1000	105	175
455 02 01	100	500	105	170
455 02 02	100	1000	105	235
455 02 03	100	2000	105	375
455 02 04	100	3000	160	300
455 02 05	100	4000	160	355
455 03 01	150	2000	160	245
455 03 02	150	3000	160	300
455 03 03	150	4000	160	355
455 03 04	150	5000	160	410
455 03 05	150	6000	180	465
455 03 06	150	10000	250	430
455 03 07	150	20000	250	650



Reaction Vessels, Cylindrical, Graduated Thermostatic Jacket with GL Connection

Catalogue Number	Flange Size (DN)	Capacity (ml)	OD/ID Diameter (mm)	Height (mm)
456 01 01	60	500	90/65	190
456 01 02	60	1000	105/80	275
456 02 01	100	500	130/105	180
456 02 02	100	1000	130/105	245
456 02 03	100	2000	130/105	375
456 02 04	100	3000	190/160	325
456 02 05	100	4000	190/160	380
456 03 01	150	2000	190/160	265
456 03 02	150	3000	190/160	320
456 03 03	150	4000	190/160	375
456 03 04	150	5000	190/160	430
456 03 05	150	6000	190/160	485

