

Reaction Vessels, Cylindrical, Thermostatic Jacket Graduated withdrawal Valve and GL Connection

Catalogue Number	Flange Size (DN)	Capacity (ml)	OD/ID Jacket Diameter (mm)	Height (mm)	Connector (GL)	Valve (mm)
457 01 01	60	500	90/65	355	14	10
457 01 02	60	1000	105/80	350	14	10
457 02 01	100	500	130/105	275	18	10
457 02 02	100	1000	130/105	340	18	10
457 02 03	100	2000	130/105	470	18	10
457 02 04	100	3000	190/160	420	18	10
457 02 05	100	4000	190/160	475	18	10
457 03 01	150	2000	190/160	370	18	10
457 03 02	150	3000	190/160	415	18	10
457 03 03	150	4000	190/160	470	18	10
457 03 04	150	5000	190/160	525	18	10
457 03 05	150	6000	190/160	580	18	10



Reaction Vessels Lid, Single Neck

Catalogue Number	Flange Size (DN)	Neck Size
461 01 01	60	14/23
461 01 02	60	24/29
461 01 03	60	29/32
461 02 01	100	24/29
461 02 02	100	29/32
461 02 03	100	34/35
461 02 04	100	45/40
461 04 01	150	29/32
461 04 02	150	34/35
461 04 03	150	45/40



Reaction Vessels Lid, Two Necks, Centre Neck and Oblique Side Neck

Catalogue Number	Flange Size (DN)	Centre Neck Size (NS)	Side Neck Size (NS)
462 01 01	60	14/23	14/23
462 01 02	60	24/29	14/23
462 01 03	60	29/32	14/23
462 02 01	100	24/29	14/23
462 02 02	100	24/29	19/26
462 02 03	100	29/32	14/23
462 02 04	100	29/32	19/26
462 02 05	100	29/32	24/29
462 02 06	100	29/32	29/32
462 02 07	100	34/35	24/29
462 02 08	100	45/40	29/29
462 04 01	150	29/32	24/29
462 04 02	150	29/32	29/32
462 04 03	150	45/40	29/32

