



Reaction Vessels Lid, Three Necks, Centre Neck and Two Oblique Side Neck

Catalogue Number	Flange Size (DN)	Centre Neck Size (NS)	Side Neck Size (NS)
463 01 01	60	14/23	14/23
463 01 02	60	24/29	14/23
463 01 03	60	29/32	14/23
463 02 01	100	24/29	14/23
463 02 02	100	24/29	19/26
463 02 03	100	29/32	14/23
463 02 04	100	29/32	19/26
463 02 05	100	29/32	24/29
463 02 06	100	29/32	29/32
463 02 07	100	45/40	29/32
463 04 01	150	29/32	24/29
463 04 02	150	29/32	29/32
463 04 03	150	45/40	29/32



Reaction Vessels Lid, Four Necks Centre Neck, Two Oblique Side Necks and One Vertical Side Neck

Catalogue Number	Flange Size (DN)	Centre Neck Size	Oblique Neck Size	Vertical Neck Size
464 01 01	60	14/23	14/23	14/23
464 01 02	60	29/32	14/23	14/23
464 01 03	60	29/32	14/23	14/23
464 02 01	100	29/32	24/29	14/23
464 02 02	100	29/32	29/32	14/23
464 02 03	100	45/40	29/32	14/23
464 04 01	150	29/32	29/32	14/23
464 04 02	150	45/40	29/32	14/23



Reaction Vessels Lid, Flat Flange Five Necks Centre Neck, Three Oblique Side Necks and One Vertical Size Neck

Catalogue Number	Flange Size (DN)	Centre Neck Size	Oblique Neck Size	Vertical Neck Size
465 02 01	100	29/32	14/23	14/23
465 02 02	100	29/32	29/32	14/23
465 02 03	100	45/40	29/32	14/23
465 04 01	150	29/32	29/32	14/23
465 04 02	150	45/40	29/32	14/23

